

PUSAT PENGURUSAN MAKMAL UNIVERSITI (PPMU)

Form Num.	UIRL/F/16
Version	2/2025
Effective Date	15/05/2025
Equipment	LPSA
Sample Serial No.	UIRL/
Page	1 of 3

PARTICLE SIZE ANALYZER LABORATORY

SAMPLE SUBMISSION FORM

General Rules and Requirements:

G CC.	ii itaic.	s and Requirements :						
1.	All information provided should be true.							
2.	Sample submission procedure.							
	a. Complete the Sample Submission Form including a valid research VOT number.							
	b. For sample submission via walk-in: Submit the completed Sample Submission Form and samples to Sample Acceptance Counter							
	c. For sample submission via mail: Submit the completed Sample Submission Form and the samples. Samust be packaged in a suitable container for courier delivery. The parcel should be addressed to the point charge of the instrument, as it will be received directly by them.							
3.	Fast lane Service: A priority testing service that provides results within 3 to 7 working days instead of the usual 14 working days. It is offered based on availability with an additional 50% charge from the normal price. Customers must contact the person in charge for this service.							
4.	For sample criteria and conditions, refer to UIRL Sample Submission Criteria in the PPMU website at ppmu.utm.my.							
5.	PPMU has the right to cancel any analysis if the sample is suspected to have a high risk on the safety of the operator or can cause damage to the instrument during the analysis. The cost of damages will be borne by the customer. Posted samples will be received by laboratory personnel.							
6.	Only samples that are ready to be analyzed are accepted by the lab.							
7.	The remaining samples will be disposed of within a month after analysis is completed.							
8.	Quotation will be provided upon request.							
9.	Payment must be made within fourteen (14) working days after the invoice is issued.							
10.	Analysis duration is within fourteen (14) working days after receiving the samples.							
11.	The laboratory will provide test results after the payment proof is presented to the laboratory personnel.							
12.	All inquiries regarding the Laser Particle Size Analyzer should be forwarded to Ms. Athirah Hanis Maulat Dzulkapli (email: athirah@utm.my) or Mrs. Shamimi Ismail (email: shamimi.ismail@utm.my) or visit our website at ppmu.utm.my							

^{*}All pages must be submitted



PUSAT PENGURUSAN MAKMAL UNIVERSITI (PPMU)

Form Num.	UIRL/F/16
Version	2/2025
Effective Date	15/05/2025
Equipment	LPSA
Sample Serial No.	UIRL/
Page	2 of 3

PARTICLE SIZE ANALYZER LABORATORY

SAMPLE SUBMISSION FORM

Application Details:

1. APPLICANT'S PERSONAL PARTICULARS								
Name of Applicant								
Status of Applicant	Und	lergraduates	Master		PhD			Researcher
Student Matric No.							•	•
Faculty/ Department								
Phone No.								
Email								
2. SUPERVISOR DETAILS (for internal applicant and academic institution only)								
Name of Supervisor	lame of Supervisor							
Staff ID No.								
Faculty/Department								
Phone No.								
Email								
	*A digital signature is not recommended. Any matters raised in the future are beyond our responsibilitie							
Signature & Official Stamp								
	I have read and agreed to the General Rules and Requirements							
3. PAYMENT								
Method of Payment	UTM PayHub System			Log card		Invoic	e	
Mode of Service		Normal Fast Lane						
Payment using Invoice	Research Vot No. (e.g.: Q.J091600.24C3.01D32)							
.,	Balance of V29000							



PUSAT PENGURUSAN MAKMAL UNIVERSITI (PPMU)

Form Num.	UIRL/F/16		
Version	2/2025		
Effective Date	15/05/2025		
Equipment	LPSA		
Sample Serial No.	UIRL/		
Page	3 of 3		

PARTICLE SIZE ANALYZER LABORATORY

SAMPLE SUBMISSION FORM

4. SAMPLE INFORMATION									
Sample Label & Information									
Sample Properties (Please tick (/))		Toxic		Carcinogenic			Others:		
Sample Type		Powder		Liqu	id	·			
			Wet Dispersion			Hydro Sight (Wet Dispersion ONLY)			
Measurement Type		Dry Dispersion		Small Volume (SV)					
				Exchangable Volume (EV)					
Particle Type		Spherical		Non-Spherical			paque Particle raunhofer Approximation)		
Sample Material, Refractive Index	Material :								
& Absorption *If no data given, default material is	Refractive Index :								
Polystyrin Latex (RI : 1.590 Absorption : 0.010)	Absorption:								
Obscuration Limit (%)	Lower:			Higher:					
Ultrasonic (if needed)	Time: mins			Temperature :°C					
Expected Result									