



UTM
UNIVERSITI TEKNOLOGI MALAYSIA

**PUSAT PENGURUSAN MAKMAL
UNIVERSITI (PPMU)**

Form Num.	UURL/F/34
Version	1/2023
Effective Date	01/02/2023
Equipment	UHPLC-ELSD
Sample Serial No.	UURL/

**LIQUID CHROMATOGRAPHY LABORATORY
SAMPLE SUBMISSION FORM (INDUSTRY)**

General Rules and Requirements:

- All information provided should be true
- Booking will be notified/updated by email
- Booking procedure
 - Submit the complete application form to UURL Sample Acceptance Counter
 - Fast Lane is offered to non-UTM customers with an additional 50% charge from the normal price**
- Sample Condition & Preparation
 - PPMU has the right to cancel any analysis if the sample is suspected to have a high risk on the safety of the operator or can cause damage to the instrument during the analysis. The cost of damages will be borne by the customer.**
 - The remaining samples will be disposed of within a month after analysis is completed.**
 - Only samples that were ready to be analyzed were accepted by the lab.
- All enquiries regarding **UHPLC** should be forwarded to the Science Officer (Mr. Ahmad Muslehuddin Sarun, email: a.muslehuddin@utm.my, tel: 07-5557775 or Assistant Engineer, Mr. Amirul Amin Khir Anuar, email: amirulamin@utm.my, tel: 07-555720) or visit our website at ppmu.utm.my

1. APPLICANT'S PERSONAL PARTICULARS									
Name of Applicant									
Hand Phone No									
Email									
Department/Division									
Signature & Official Stamp		*A digital signature is not recommended. Any matters raised in the future are beyond our responsibilities							
2. COMPANY DETAILS									
Name									
Registration No.									
Address									
Telephone No.									
Email									
Mode of Payment		<input type="checkbox"/>	Cash	<input type="checkbox"/>	EFT	<input type="checkbox"/>	Invoice	<input type="checkbox"/>	Fast Lane
3. SAMPLE INFORMATION									
Total No. of Sample									
Name of Sample/s									
Sample Properties (please tick (✓))		<input type="checkbox"/>	Toxic	<input type="checkbox"/>	Carcinogenic	<input type="checkbox"/>	<input type="checkbox"/>	Normal	<input type="checkbox"/>
Sample Purity									



UTM
UNIVERSITI TEKNOLOGI MALAYSIA

**PUSAT PENGURUSAN MAKMAL
UNIVERSITI (PPMU)**

Form Num.	UIRL/F/34
Version	1/2023
Effective Date	01/02/2023
Equipment	UHPLC-ELSD
Sample Serial No.	UIRL/

**LIQUID CHROMATOGRAPHY LABORATORY
SAMPLE SUBMISSION FORM (INDUSTRY)**

Targeted Compounds	
--------------------	--

4. ANALYSIS INFORMATION (please attach the copy of referred journal)

Type of Column Available	Waters, XBridge BEH C18, 4.6 x 250mm, 5 micron		Waters, XBridge Amide BEH, 4.6 x 250mm, 5 micron				
	Waters, XBridge Phenyl BEH, 4.6 x 250mm, 5 micron		ZORBAX Eclipse Plus C18, Rapid Resolution 4.6 x 100mm 3.5 micron				
	Phenomenex Gemini® 5µm NX-C18 LC column 250 x 4.6 mm						
Elution (please tick (√))	Isocratic		Gradient				
Injection Volume (µL)							
Flow Rate (mL/min)							
Gain							
Oven Temperature (°C)							
Stoptime (min)							
Postrun (min)							
Mobile Phase or Premix (If Isocratic)	A :	%	C :	%			
	B :	%	D :	%			
Mobile Phase Timetable (If Gradient)	Time (min)	A (%)	B (%)	C (%)	D (%)	Flow (mL/min)	Max Pressure (bar)